


## • Animal Handling in *In-Vivo* Scanners




### Mouse and rat cassettes

The SkyScan *in-vivo* systems are supplied with exchangeable animal cassettes that can be used in other Bruker *in-vivo* imaging instruments such as MRI, micro-PET, micro-SPECT, bio-luminescence, bio-fluorescence, etc. to collect multimodal information. Such data are used for co-registration of functional and morphological information from the same animal.

The cassettes are equipped with face mask and tubes for anaesthetic gas. There are also connections to ECG electrodes and a temperature sensor. All tubes and contacts are combined in a single connector, which can be attached to or detached from the animal transport system by a small slider. To prevent leakage of anaesthetic gas, the connections have valves to stay closed if the animal cassette is disconnected from the animal transport system.

### All-carbon ECG electrodes


The animal cassettes include special cables and clip electrodes to detect ECG signals by a sensitive amplifier integrated in the physiological monitoring sub-system. The ECG connections use special wiring and electrodes developed by Bruker microCT which contain no metal parts. The wires and electrodes employ advanced carbon-fiber conductive parts with X-ray absorption similar to that of animal tissues for uncompromised image quality.



Shadow projections (top) and reconstructed slices (bottom) through a mouse body with attached metal electrodes (left) and carbon electrodes developed by Bruker microCT (right).



## • Integrated Physiological Monitoring



### Physiological monitoring

The physiological monitoring sub-system includes video monitoring of the animal with real-time movement detection, ECG and breathing detection, and temperature stabilization. A 5 megapixel color camera is mounted above the animal bed and equipped with white LED illumination to produce a real-time image of the animal during the scan. The software analyzes the video stream from a user-selected area of the image, which the operator can position on a part of the animal body where breathing movement is visible. These movements are converted into a movement waveform to provide timemarks for time-resolved micro-CT reconstruction.

The face mask on the animal bed is connected to an air/gas flow sensor for direct breathing detection. The ECG electrodes in the animal cassette are connected to a sensitive ECG amplifier. Both breathing and ECG signals are digitized, sent to computer and displayed as real-time profiles on-screen. An operator can select individual gain and threshold for each signal to optimize generation of timemarks. The monitoring also includes temperature stabilization by heated airflow, which maintains the scanned animal at a selected temperature, to prevent cooling of the animal under anaesthesia.

### 4D time-resolved microtomography

Physiological monitoring creates reference time-marks for time-resolved reconstruction of heart and lung dynamics, invented by Bruker microCT. Using such timemarks, multiple projection images taken at each gantry angular position are sorted post-scan into breathing or heart time bins. Such sorting creates pseudo-static sets of projections, which are reconstructed as separate datasets and produce 3D sets of results corresponding to different phases of the cardiac or respiratory cycle.

Our visualization program loads reconstructed datasets and allows scrolling in XYZ dimensions across the reconstructed volume. In the time-dimension it visualizes dynamics of heart or lung movements in sharp reconstructed images minimally affected by movement artifacts. Because all acquired data are sorted after the acquisition process, respiratory and cardiac cycles can be visualized by sorting according to timemarks from physiological monitoring without rescanning the animal.

